

OriGene Technologies Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850 UNITED STATES Phone: +1-858-888-7900 Fax: +1-858-888-7904 <u>US-info@acris-antibodies.com</u>

BP740 **OriGene EU**

Acris Antibodies GmbH Schillerstr. 5 32052 Herford GERMANY Phone: +49-5221-34606-0 Fax: +49-5221-34606-11 info@acris-antibodies.com

	Polyclonal Antibody to Algea allo-Phycocyanin
Catalog No.:	BP740
Quantity:	1 ml
Host:	Rabbit
Immunogen:	Native.
	Format: This antibody is supplied as liquid, affinity purified IgG in 50 mM borate buffered saline, pH 8.0, 0.05% sodium azide as preservative.
Applications:	Immunohistochemistry. ELISA. Western Blot. Dilution: Do not dilute the entire recondtituted solution at once. Dilute with 50 mM borate buffered saline, pH 8.0, 0.02% sodium azide and 0.15M sodium chloride. Other applications not tested. Optimal dilutions of this antibody are dependent on conditions and should be determined by the user. Other applications not tested. Optimal dilutions are dependent on conditions and should be determined by the user.
Specificity:	Reacts with allo-Phycocyanin and Phycocyanin. Does not react with Phycoerythrin.
Storage:	Store at 4°C, DO NOT FREEZE. Shelf life: on

For research and in vitro use only. Not for diagnostic or therapeutic work. Material Safety Datasheets are available at www.acris-antibodies.com or on request. Acris Antibodies is now part of the OriGene family. Learn more at www.origene.com



OG/20121030