

### OriGene Technologies Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850 UNITED STATES

Phone: +1-858-888-7900 Fax: +1-858-888-7904 US-info@acris-antibodies.com

# PA1192XC OriGene EU

#### Acris Antibodies GmbH

Schillerstr. 5 32052 Herford GERMANY

Phone: +49-5221-34606-0 Fax: +49-5221-34606-11 info@acris-antibodies.com

## **Human Menopausal Gonadotropin - Purified**

Catalog No.: PA1192XC

Quantity: 10 mg

Concentration: 1.0 mg/ml

Species: Human

**Source:** Urine, Urine of post-menopausal women.

Format: State: Sterile Filtered White lyophilized (freeze-dried) powder.

Purity: >98% Greater than 98.0% as determined by:

(a) Analysis by RP-HPLC.(b) Anion-exchange FPLC.

(c) Analysis by reducing and non-reducing SDS-PAGE Silver Stained gel.

**Buffer System:** Lyophilized from a concentrated solution in water containing no additives. **Reconstitution:** It is recommended to reconstitute the lyophilized hMG in sterile  $18M\Omega$ -cm H2O not less than  $100\mu g/ml$ , which can then be further diluted to other aqueous solutions.

**Description:** hMG is produced from a sterile preparation of placental glucoprotein urine of post-

menopausal women. The hMG is purified by proprietary chromatographic techniques. **Biological Activity:** hMG conforms the following activities: 100IU/mg of FSH and 100 IU/mg

LH.

Add. Information: Free of: HbsAg, Hepatitis B surface antigen and antibodies to HIV, Hepatitis C and HCV.

Protein quantitation was carried out by two independent methods:

1. UV spectroscopy at 280 nm.

2. Analysis by RP-HPLC, using a standard solution of hMG as a Reference Standard.

Storage: Lyophilized hMG although stable at room temperature for 3 weeks, should be stored

desiccated below -18 C. Upon reconstitution hMG should be stored at 4 C between 2-7 days and for future use below -18 C. For long term storage it is recommended to add a carrier

protein (0.1% HSA or BSA). Please avoid freeze-thaw cycles.

General Readings: 1. Lin YH, Seow KM, Chen HJ, Huang LW, Hwang JL, Tzeng CR. Impact of estradiol patterns in

clomiphene citrate/human menopausal gonadotropin/cetrorelix protocol. Gynecol

Endocrinol. 2007 Jan; 23(1): 45-9. PubMed PMID: 17484512.

2. Zhang HQ, Yan B, Zhao HX, Gu DY, Jia XF, Cao L, et al. Effect of traditional Chinese herbs combined with low dose human menopausal gonadotropin applied in frozen-thawed embryo transfer. Chin J Integr Med. 2006 Dec;12(4):244-9. PubMed PMID: 17361518.

3. Lin YH, Hwang JL, Seow KM, Huang LW, Hsieh BC, Tzeng CR. Comparison of outcome of clomiphene citrate/human menopausal gonadotropin/cetrorelix protocol and buserelin long protocol--a randomized study. Gynecol Endocrinol. 2006 Jun;22(6):297-302. PubMed

PMID: 16785154.

4. Zech NH, Netzbandt P, Zech H. Development of severe ovarian hyperstimulation syndrome after inadvertent stimulation with a gonadotropin-releasing hormone agonist

TUV HORD CERT CRIENT OF THE PARK AND THE THE P



#### PA1192XC: Human Menopausal Gonadotropin - Purified

and human menopausal gonadotropin in a pre-existing early pregnancy. Fertil Steril. 2005 Dec;84(6):1745. PubMed PMID: 16359980.

- 5. Keye WR, Webster B, Dickey R, Somkuti S, Crain J, Scobey MJ. Subcutaneously administered Menopur, a new highly purified human menopausal gonadotropin, causes significantly fewer injection site reactions than Repronex in subjects undergoing in vitro fertilization. Reprod Biol Endocrinol. 2005 Nov 9;3:62. PubMed PMID: 16280073.
- 6. Ov ulation induction in the new millennium: recombinant follicle-stimulating hormone versus human menopausal gonadotropin. Gynecol Endocrinol 2005 Mar;20(3):161-9