

NEUROMICS



Neurofibrillary Tangles

Data Sheet

Catalog Number: RA19077

Host: Rabbit

Product Type: Whole Serum Antibody

Species: Human

Reactivity:

Immunogen Sequence: Disaggregated paired helical filaments were used as an immunogen. Insoluble NFT which were extracted in boiling sodium dodecyl sulfate, purified by sucrose gradient centrifugation and disaggregated by sonication were used as an immunogen for the preparation of this antibody.

Format: Liquid. Whole Serum

Applications: Western blot: 1:500
Immunohistochemistry: 1:200 (formalin-fixed, paraffin-embedded sections)
1:2,000 (paraformaldehyde-fixed, frozen sections).

Dilutions listed as a recommendation. Optimal dilution should be determined by investigator.

Storage: Maintain at +2-8°C for 3 months or at -20°C for longer periods. Stable for 1 year. *Avoid repeated freeze-thaw cycles.*

References: Skrabana, R., Skrabanova-Khuebachova M., Kontsek, P. and Novak M. Alzheimer's-disease-associated conformation of intrinsically disordered tau protein studied by intrinsically disordered protein liquid-phase competitive enzyme-linked immunosorbent assay. *Anal Biochem.* 2006 Dec 15;359(2):230-7.

Arima, K. Ultrastructural characteristics of tau filaments in tauopathies: immuno-electron microscopic demonstration of tau filaments in tauopathies. *Neuropathology.* 2006 Oct;26(5):475-83.

Application Notes

This antibody specifically immunostains NFT and degenerating plaque neurites in cases of AD, Down's syndrome and normal aged individuals. NFT are heterogenous in composition and their composition undergoes changes as the tangle matures. This antibody labels all NFT, both intracellular and extracellular, as well as all paired helical filaments (PHF) containing neurites.

FOR RESEARCH USE ONLY

NEUROMICS REAGENTS ARE FOR IN VITRO AND CERTAIN NON-HUMAN IN VIVO EXPERIMENTAL USE ONLY AND NOT INTENDED FOR USE IN ANY HUMAN CLINICAL INVESTIGATION, DIAGNOSIS, PROGNOSIS, OR TREATMENT. THE ABOVE ANALYSES ARE MERELY TYPICAL GUIDES. THEY ARE NOT TO BE CONSTRUED AS BEING SPECIFICATIONS. ALL OF THE ABOVE INFORMATION IS, TO THE BEST OF OUR KNOWLEDGE, TRUE AND ACCURATE. HOWEVER, SINCE THE CONDITIONS OF USE ARE BEYOND OUR CONTROL, ALL RECOMMENDATIONS OR SUGGESTIONS ARE MADE WITHOUT GUARANTEE, EXPRESS OR IMPLIED, ON OUR PART. WE DISCLAIM ALL LIABILITY IN CONNECTION WITH THE USE OF THE INFORMATION CONTAINED HEREIN OR OTHERWISE, AND ALL SUCH RISKS ARE ASSUMED BY THE USER. WE FURTHER EXPRESSLY DISCLAIM ALL WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. v1-04/10

www.neuromics.com

Neuromics Antibodies • 5325 West 74th Street, Suite 8 • Edina, MN 55439
phone 866-350-1500 • fax 612-677-3976 • e-mail pshuster@neuromics1.com