

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850 UNITED STATES Phone: +1-888-267-4436 Fax: +1-301-340-8606 techsupport@origene.com **OriGene Technologies GmbH** 

Schillerstr. 5 32052 Herford GERMANY Phone: +49-5221-34606-0 Fax: +49-5221-34606-11 info-de@origene.com

## SM3164P

**Quantity:** 

Alternate names:

**Concentration:** 

**Background:** 

## Monoclonal Antibody to alpha + beta Tubulin (heterodimer) - Ig Fraction

alpha-tubulin, beta-tubulin, tubulin alpha, tubulin beta 0.1 mg 1.0 mg/ml The microtubules are intracellular dynamic polymers made up of evolutionarily conserved polymorphic alpha/beta-tubulin heterodimers and a large number of microtubule-associated proteins (MAPs). The microtubules consist of 13 protofilaments and have an outer diameter 25 nm. Microtubules have their intrinsic polarity; highly dynamic plus ends and less dynamic minus ends. Microtubules are required for vital processes in eukaryotic cells including mitosis, meiosis, maintenance of cell shape and intracellular transport. Microtubules are also necessary for movement of cells by means of flagella and cilia. In mammalian tissue culture cells microtubules have their minus ends anchored in microtubule organizing centers (MTOCs). The GTP (guanosintriphosphate) molecule is an essential for tubulin heterodimer to associate with other heterodimers to form microtubule. In vivo, microtubule dynamics vary considerably. Microtubule polymerization is reversible and a populations of microtubules in cells are on their minus ends either growing or shortening - this phenomenon is called dynamic instability of microtubules. On a practical level, microtubules can easily be stabilized by the addition of nonhydrolysable analogues of GTP (eg. GMPPCP) or more commonly by anti-cancer drugs such as Taxol. Taxol stabilizes microtubules at room temperature for many hours. Using limited proteolysis by enzymes both tubulin subunits can be divided into Nterminal and C-terminal structural domains. The alpha-tubulin (relative molecular weight about 50 kDa) is globular protein that exists in cells as part of soluble alpha/beta-tubulin dimer or it is polymerized into microtubules. In different species it is coded by multiple tubulin genes that form tubulin classes (in human 6 genes). Expressed tubulin genes are named tubulin isotypes. Some of the tubulin isotypes are expressed ubiquitously, while some have more restricted tissue expression. Alphatubulin is also subject of numerous post-translational modifications. Tubulin isotypes and their posttranslational modifications are responsible for multiple tubulin charge variants - tubulin isoforms. Heterogeneity of alpha-tubulin is concentrated in Cterminal structural domain. The beta-tubulin (relative molecular weight about 50 kDa) is counterpart of alpha-tubulin in tubulin heterodimer, it is coded by multiple tubulin genes and it is also posttranslationally modified. Heterogeneity of subunit is concentrated in C-terminal structural domain. The beta-tubulin (relative molecular weight about 50 kDa) is counterpart of alphatubulin in tubulin heterodimer, it is coded by multiple tubulin genes and it is also posttranslationally modified. Heterogeneity of subunit is concentrated in C-terminal structural domain.

Host / Isotype:

Mouse / IgM

For research and in vitro use only. Not for diagnostic or therapeutic work. Material Safety Datasheets are available at www.acris-antibodies.com or on request.

OG/20130815 **1** 

SM3164P: Monoclonal Antibody to alpha + beta Tubulin (heterodimer) - Ig
Fraction

	Fraction
Clone:	TU-08
Immunogen:	Microtubule proteins from porcine brain
Format:	<b>State:</b> Liquid Ig fraction <b>Purification:</b> Precipitation methods and size-exclusion chromatography <b>Buffer System:</b> PBS with 15 mM sodium azide, approx. pH 7.4
Applications:	Western blotting: 5 mg/ml. Positive control: porcine brain lysate. Sample preparation: Mix lysate with reducing Laemmli SDS-PAGE sample buffer. Application note: Reducing conditions. Immunocytochemistry. Other applications not tested. Optimal dilutions are dependent on conditions and should be determined by the user.
Specificity:	The antibody recognizes alpha, beta-tubulin heterodimer (porcine brain), a basic intracellular structural unit of microtubules. Alpha- (50 kDa) and beta- (50 kDa) tubulins form tubulin heterodimer, a globular protein that polymerizes to form microtubules. <b>Species:</b> Human, Porcine, Mouse. Other species not tested.
Storage:	Store at 2-8°C. Do not freeze. Shelf life: one year from despatch.
General Readings:	<ol> <li>Linhartova I, Draberova E, Viklicky V, Draber P.: Distribution of non-class-III beta- tubulin isoforms in neuronal and non-neuronal cells. FEBS Lett. 1993 Mar 29;320(1):79-82.</li> <li>Cowan NJ, Dobner PR, Fuchs EV, Cleveland D.W.: Expression of human alpha- tubulin genes: interspecies conservation of 3' untranslated regions.Mol Cell Biol. 1983 Oct;3(10):1738-45.</li> <li>Draberova E, Draber P.: Novel monoclonal antibodies TU-08 and TU-16 specific for tubulin subunits. Folia Biol (Praha). 1998;44(1):35-6.</li> </ol>

For research and in vitro use only. Not for diagnostic or therapeutic work. Material Safety Datasheets are available at www.acris-antibodies.com or on request.